

# Active oxygen granules

## For gentle chlorine-free water disinfection

#### **Application**

Active oxygen granules will disinfect your spa water gently. This totally chlorine-free treatment method creates a more natural and pleasant water quality. To increase the efficacy, SpaTime active oxygen activator must be used along with active oxygen granules.

**Advantages** 

- Pleasant and odourless water quality without chlorine
- Gentle to skin and hair
- Free of any metals
- Fast dissolving
- Dissolves completely without residue

### **Product description**

Fast dissolving, white active oxygen granular Contains: 100 g of the product contains 100 g Pentapotassium bis(peroxymonosulphate) bis(sulphate)

### **Dosing recommendations**

Check the active oxygen content with SpaTime test strips as soon as the granules have completely dissolved, as it will not be possible to obtain a correct measurement at a later time. The content should be at least 8 mg/L. If it is not, add the same amount again. Before each use, add the active oxygen granules directly to

the spa water with the circulation pump running. If the spa is used frequently, you need to double the amount in order to disinfect it sufficiently. The active oxygen granules only reach optimum effectiveness when used along with the active oxygen activator.

#### **Important information**

This product should be used solely for the purpose specified.

TIP: Before adding the disinfectant, please check the pH value with SpaTime test strips. If necessary, adjust to the correct value of 7.0 – 7.6 to ensure the best possible efficacy of the active oxygen granules.



